

Only \$340

Usually \$484

Code: CORN6770

A personal benchtop instrument designed for quick spin downs of micro-samples.



- Quick release rotor system eliminates the need for tools to exchange rotors
- Electronic braking provides rapid deceleration for more efficient sample prep
- Compact, rugged design withstands heavy use in crowded lab environment
- Dual lid switch provides safe operation
- Ergonomic design allows for comfortable handling
- Includes 8 x 1.5/2.0mL rotor and 4-place PCR strip tube rotors
- Fixed speed of 6,000rpm generates 2,000 x g