

## Positive Displacement Dispenser Tips

**Buy 2 Boxes Get  
7% Off**

100 tips per box

An efficient tool for quick and precise filling of plates or large series of tubes as used in nucleic acid purification kits.

Optimize your workflows involving plates or large series of other vessels with Combitips advanced dispenser tips and accurately and precisely pipette both viscous and volatile liquids.

These positive displacement pipette tips are your consumable of choice for all Multipette® multi-dispenser pipettes from Eppendorf and offer a variety of advantages in comparison to their predecessor Combitips plus.

Available in nine volume sizes (0.1 mL – 50 mL), Combitips advanced have a maximum range of dispensing volumes with increments as low as 100 nanoliters depending on the Multipette model used.

- Available in 4 purity grades including Eppendorf Sterile, Quality, Biopur® and Forensic DNA Grade
- The elongated tips of the Combitips advanced make it possible to empty all standard lab tubes entirely
- Individual color-coding makes it easy to quickly identify the desired Combitips advanced



See page 2 for available tips

# Combitips® Advanced

Code	Size	Colour
------	------	--------

## Eppendorf Quality

EPPE0030089.413	0.2ml	Light Blue
EPPE0030089.421	0.5ml	Purple
EPPE0030089.464	10ml	Orange
EPPE0030089.480	50ml	Light Grey

## Biopur®

EPPE0030089.618	0.1ml	White
EPPE0030089.626	0.2ml	Light Blue
EPPE0030089.634	0.5ml	Purple
EPPE0030089.650	2.5ml	Green
EPPE0030089.669	5ml	Blue
EPPE0030089.677	10ml	Orange
EPPE0030089.685	25ml	Red
EPPE0030089.693	50ml	Light Grey

## PCR Clean

EPPE0030089.782	0.5ml	Purple
-----------------	-------	--------

