



# Stirrer Magnetic C-MAG HS 7 Digital Package



**Only \$1,695**

Usually \$2,260

Code:IKAA10010600

### Package Includes:

- |                        |                          |
|------------------------|--------------------------|
| 1 x C-MAG HS 7 digital | 1 x H 38 Holding rod     |
| 1 x H 16 V Support rod | 1 x H 44 Boss Head Clamp |

Magnetic stirrer with heating and ceramic heating plate which offers excellent chemical resistance

- Powerful motor for stirring quantities of up to 10 l (H<sub>2</sub>O)
- Simultaneous digital display of target and actual temperatures via LCD display
- Directly connection for a PT 1000 temperature sensor enable a precise temperature control (included in delivery)
- Control accuracy in medium +/- 0,5 K (in combination with PT 1000)
- Display of actual temperature value in medium with 0,1 K resolution when using PT 1000 temperature sensor
- 3 operating modes to choose (standard, safe, adjustment protection)
- Fixed safety circuit of 550 °C
- Hot Top indicator >> hot surface warning to prevent burns!
- Digital error code display
- Elevated control panel for protection against leaking liquids



Ref. Code:  
IKA-HS7-PK

VALID UNTIL:  
31 July 2024



1800 00 84 53

[sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

[shop.bio-strategy.com](http://shop.bio-strategy.com)