



# Miele

## Professional PG8504 Washer



**Only \$7,833**

Usually \$8,704

Code: MIEL62850422

Remove the monotony & hazard potential of washing glassware manually, let's put some precision and safety into this procedure!

By removing the manual washing step from your wash process there is reduced likelihood of glass breakage that can result in cuts to lab staff.

- Powerful 500L/min pump with integrated heating elements enabling safe washing cycles with temperatures up to 70°C
- Underbench/Free-Standing
- Door dispensers for powder detergent and liquid neutraliser
- Suitable for operation on 240V 50Hz 13 amp supply
- 6 standard programs + 2 vacant program positions to meet your wash requirements
- Touch on steel controller with 3-line display for easy viewing of washer status
- Compact: 600mm x 600mm x 850mm W x D x H
- Stainless steel interior/exterior for durability and cleanliness



Miele - good for occupational health and safety

Many options and accessories to select from.

Contact Bio-Strategy for assistance in configuring your washer.



**bio-strategy**  
Part of DKSH Group

Ref. Code:  
MIEL-PG8504

VALID UNTIL:  
31 Aug 2024



**DKSH**

 1800 00 84 53

 [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)

 [www.bio-strategy.com](http://www.bio-strategy.com)

 [shop.bio-strategy.com](http://shop.bio-strategy.com)

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.