



725L FrostLess -86°C ULT Freezer

Combat frost build-up damage to
freezer components with this low frost model



Introductory Offer!

Contact us for a quote today

Code: SANYMDF-DU700ZH-PE

Temperature Range: -40°C to -86°C

The FrostLess ultra-low freezer ensures the reliability and energy efficiency modern laboratories demand with less frost build-up.

The cutting-edge thermal insulation design helps combat frost build-up to reduce freezer stress and preserve high performance.

Meet Sustainability and Budget Goals

The use of energy efficient, natural refrigerants combined with an inverter-controlled compressor, lowers operational costs, environmental impact, and total cost of ownership.



Ref. Code:
SANY-Frostless

VALID UNTIL:
30 Sep 2024



1800 00 84 53

sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.