

# Incubate, Shake, Or Do Both

## INC 125 FS Digital Incubator Shaker

The first incubator shaker with a patented removable shaker platform



**watch >**

**Only \$9,495**

Usually \$12,665

Code:	Orbit
IKAA20117226	20mm
IKAA20117235	25mm

Use as an incubating shaker with shaker platform, or as an incubator using the wire rack (included).

- High shaking speed up to 300 rpm
- Orbit: 20 mm and 25 mm models available
- Wide temperature range: Ambient +8°C to 80 °C
- Adjustable fan speed (0-100%)
- Corrosion-resistant stainless-steel interior
- Decontamination mode for disinfection and reduction of cross-contamination (100/120°C)
- Counter, timer and timer auto (temperature-controlled) function
- 2-year warranty



Space-saving stackable design

Additional accessories available (flask clips, sticky mats, etc) to accommodate a wide range of vessels



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
IKAA-INC125FS

VALID UNTIL:  
31 July 2024