

## PG 8583 CD ADP Washer/Disinfector

**Only \$19,580**  
Usually \$21,755  
Code: MIEL62858351



Lab washer with pump for AD water, liquid dispensing, DryPlus and conductivity monitoring.  
Have peace of mind with German quality and a 2-year warranty.

- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- 14 standard programs + 2 freely configurable programs
- Two wash levels configured for washing wide neck and narrow neck glassware
- Conductivity monitoring to ensure load is clean
- Dispensing pumps for liquid detergent and liquid neutraliser
- Spray arm and spray pressure monitoring ensuring a clean wash
- Underbench or standalone model
- Stainless steel interior/exterior for durability and cleanliness
- ADP model including pump for AD water and waterproofing kit
- Many options and accessories to select from
- Pull out drawer for housing detergent and neutraliser

Complete system requires optional baskets and inserts.

Please contact Bio-Strategy for assistance in configuring your complete washing solution.

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
MIEL-PG8583CD

VALID UNTIL:  
31 July 2024