

AZLON® Ethanol Wide Neck Wash Bottle

AZLON® multi-lingual safety wash bottles ensure the safe and accurate dispensing of liquids.

Available in 11 languages and with 6 solvent labels, you can be sure to find the perfect wash bottle for you.



NOTE:

For dispensing use only. Not for liquid storage.

Other wide neck bottles also available

Vented

With the patented Driplok® vapour venting valve that prevents hazardous solvent drips.

Code	Usually	Now
	5 bottles per pack	
AZLOWGW535VTML 250ml	\$55.80	\$44.60
AZLOWGW541VTML 500ml	\$58.40	\$46.70

Non-Vented

Flexible LDPE bottle fitted with a single piece colour coded polypropylene closure

Code	Usually	Now
	5 bottles per pack	
AZLOWGW541PML 500ml	\$47.70	\$38.20

Buying in bulk?

Contact us and we can strike a deal!

- Manufactured from flexible LDPE
- Indelibly printed in up to four languages per bottle
- Closure features precision jet-tip aperture for fine control
- Includes the new Globally Harmonised Symbols, NFPA diamond, risk phrases and safety information