

Pyrex Tall Form Beaker with Spout

**SAVE
20%**

Code	Size	Was	Now
BIBB1015/02M	50ml	\$53.10	\$42.50
BIBB1015/06D	150ml	\$63.10	\$50.50
BIBB1015/10D	250ml	\$66.20	\$53.00
BIBB1015/14D	400ml	\$71.60	\$57.30

Buying in bulk?

Contact us and we can strike a deal!

Suitable for use in a range of laboratory applications including titrations.

Pyrex® beakers are printed with Trace Code allowing a batch certificate to be downloaded.

- Complies with ISO 3819 and DIN 12331
- Manufactured from borosilicate glass 3.3
- White graduations and marking spot
- Graduations are approximate

