

ScanSpeed MiniVac Beta

Simple, Straight-Forward & Trouble-Free
Speed Vacuum Concentration

Only \$13,400

Usually \$16,335

Code: LABO7008100778

Designed for rapid and efficient concentration
or drying of small samples in tubes.

Supplied complete with Scanspeed 32 concentrator,
multipurpose rotor, Vacuubrand membrane pump,
DC drip catcher, ice cube container, and all connectors,
hoses, and clips.



Digital display of the vacuum down to 1 mBar, full programming of time and temperature.

The performance can be enhanced by addition of ice cubes around the drip catcher to increase the rate of evaporation.

Acrylic lid allows full visibility of the ScanSpeed 32 centrifuge chamber to monitor the progress.

No cross contamination or sample loss, the exhaust vent for the vacuum being totally isolated.

- Capacity: 48 x 1.5-2.0 mL / 76 x 0.5mL tubes
- Speed: 200-2000 RPM
- Chamber temperature setting: Off, 4°C to 65°C
- Timer with Start/Stop
- Digital display of vacuum
- Programming of time and temperature