

Pocket-Chek Refractometers

**Compact, Powerful,
Dependable, Accurate**

Powerful multi-scale performance in
the palm of your hand!

Only \$731

Usually \$860

BRIX/RI - Code REIC13940000

Three measurement scales are offered in this unit: Refractive Index and Brix with automatic temperature compensation, and Refractive Index at actual temperature.

SALT - Code: REIC13940020

Suitable for numerous industries including Food, Beverage, Agriculture, and Industrial. Three measurement scales in one model allow more measurement flexibility and cost savings for you.

VET360 - Code: REIC13940018

Specifically designed for the veterinary market with measuring scales based on animal chemistries. Four measurement scales are provided: Canine, Feline and large animal urine specific gravity, and Serum protein concentration for all animals.

Multi - Code REIC13940015

This 4 in 1 tool features four measuring scales: Ethylene Glycol Propylene Glycol and Windshield wash freeze point, Battery electrolyte specific gravity.

GLYCERIN (Celsius) - Code: REIC13940023

This unit features three measuring scales: Freeze Point, % concentration and Boiling Point of Glycerine.

