


KIMBLE® Sterile Disposable Glass Pasteur Pipets



Only \$125

Usually \$179

Code: KIML63A53

146mm | 5.75in
1000 items

Only \$175

Usually \$219

Code: KIML63B92

228mm | 9in
1000 items

These unplugged pasteur-type pipets are intended for
one-time use in the transfer of small volumes

- Designed to ASTM E732 Standard Specification for Disposable Pasteur-Type Pipette
- Ideal for general liquid transfer applications in the laboratory
- Borosilicate pipets are manufactured from 51 expansion, low extractable borosilicate glass tubing conforming to USP Type 1 and ASTM E438, Type 1, Class B requirements
- Supplied sterile by gamma-irradiation
- Also available cotton plugged