

VACSTAR Control

Just \$3,022

Usually \$6,044

Code: IKAA20019070

1 Only



The VACSTAR dry 4-chamber membrane vacuum pump is powerful, resistant to chemicals and low-maintenance.

An ideal system with its integrated vacuum controller and removable display. Offering a suitable mode for the desired process: from automatic to manual operation, right through to pumping or an evacuation and cleaning function

Integrated vacuum controller with removable controller (WiCo)

Depending on the application, you can easily select your operating mode either by using the intuitive operator menu or by remote control.

Color Squid White

Just \$495

Usually \$785

Code: IKAA3671000

Made from recycled materials



Small magnetic stirrer with glass surface.

- Digital speed display (LED)
- Electronically controlled motor for more capacity
- Higher speed range from 0 - 2.500 rpm
- Max. stirring quantity 1 l
- Outstanding chemical resistance due to glass top and synthetic bottom made of TPC-ET

Don't miss out - Order Today!