

Introducing the VIP ECO -86°C Chest Freezers

180L
Code: SANYMDF-DC202VH-PE



[watch >](#)

Providing maximum sample storage capacity within an optimum footprint combined with industry leading cooling performance and reliability for an ideal long-term preservation solution.

Optimum Sample Storage

The use of patented VIP PLUS panels within the freezer cabinet provides a reduced wall thickness for maximum interior volume. This technology allows the unit to require a small installation space with a great storage capacity.

Minimal Footprint

Up to 30% more storage capacity than a conventionally insulated freezer without increasing the footprint with VIP PLUS vacuum panels for increased energy efficiency and cooling performance.

Give us a call for a great introductory price

84L
Code: SANYMDF-DC102VH-PE

