



# Introducing the TwinGuard ECO Pharmaceutical Refrigerator



With sliding doors it perfectly preserves samples, reagents, and pharmaceuticals through high-reliability temperature control.

### Ultimate Sample Protection

A single unit incorporates two independent refrigerating circuits. Should one circuit experience a problem, the other circuit will continue to maintain in-chamber temperature at  $SV 5^{\circ}C \pm 3^{\circ}C$ . The double-safe circuits ensure preservation of crucial stored items with greater protection.

### Improved Sustainability

Natural HC refrigerants, inverter-controlled compressors, new Cycle Defrost system, together with enhanced sliding glass doors design, combine to improve temperature management, increase energy efficiency and minimal footprint.

### User-friendly Design

The ergonomic design provides a clear view of stored items through the large-framed windows. The slim profile allows easy retrieval, and the self-closing sliding doors makes them ideal for limited spaces.

Code: SANYMPR-S1201XH-PE

Effective capacity 1165L

Temperature range 2°C to 14°C

Includes 12 hard steel wire shelves with polyethylene coating

Also available in a ½ drawer shelving format

**Give us a call for a great introductory price**



Part of DKSH Group

Ref. Code:  
SANY-ECOTwin

VALID UNTIL:  
31 Oct 2024



1800 00 84 53

sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com