

# Pyrex® Berzelius Beaker Tall Form with Spout

Also known as Berzelius beakers, tall form beakers are suitable for use in a range of laboratory applications including titrations. Pyrex® beakers are printed with Trace Code allowing a batch certificate to be downloaded.



- Compliant with ISO 3819 and DIN 12331
- Manufactured from borosilicate glass 3.3
- Printed with Trace Code for downloadable batch certificate
- White graduations and marking spot
- Graduations are approximate

BIBB1015/02M	Beaker Tall Form, 50mL	10 items	\$ 53.10
BIBB1015/06D	Beaker Tall Form, 150mL	10 items	\$ 63.10
BIBB1015/10D	Beaker Tall Form, 250mL	10 items	\$ 66.20
BIBB1015/14D	Beaker Tall Form, 400mL	10 items	\$ 71.60

**In Stock!**

**Special \$ 40**

**Special \$ 47**

**Special \$ 50**

**Special \$ 54**