



# Digital Magnetic Stirrer - RO 5 - 1 only!



**RO-5 - priced to sell**

**Only \$ 1,695**

Usually \$ 3,390

Code: IKAA3690500

5-position digital magnetic stirrer without heating, designed for synchronous stirring. The magnetic coil technology provides noiseless and consistent stirring on all positions. Speed remains constant, even when load changes.

The closed and compact design allows easy cleaning and protects the equipment against the penetration of liquids. Additionally, a removable transparent anti-slip mat is included.

- Adjustable reverse operation
- Wear free magnetic coils
- Error code display
- Easy operation with touch keypad
- Suitable for continuous operation
- Speed range: 0 - 1200 rpm



Part of DKSH Group

Ref. Code:  
IKAA\_RO5

VALID UNTIL:  
31.10.24  
current stock only



1800 00 84 53

sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



# Digital Magnetic Hotplate Stirrer RT 5 - 1 only!



**RT-5 - priced to sell**

**Only \$ 3,030**

Usually \$ 6,060

Code: IKAA3690605

5-position digital magnetic hotplate stirrer, designed for synchronous heating and stirring.

The magnetic coil technology provides noiseless and consistent stirring on all positions. The surface temperature is infinitely adjustable up to 120 °C, producing a maximum medium temperature of 70 °C (depending on the used vessel). Speed remains constant, even when load changes.

- Homogeneous temperature distribution across the heating plate
- Adjustable reverse operation
- Eco-mode for a low self-warming of the surface
- Wear free magnetic coils
- Error code display
- Easy operation with touch keypad
- Suitable for continuous operation
- Speed range: 0 - 1000 rpm
- Silicone mat for protection and optimized temperature distribution included



Part of DKSH Group

Ref. Code:  
IKAA\_RO5

VALID UNTIL:  
31.10.24  
current stock only



1800 00 84 53



sales.au@bio-strategy.com



www.bio-strategy.com



shop.bio-strategy.com