



# New and Improved IKA RCT Basic


**In Stock and  
Ready to Ship**

**Only \$1,195**

**Usually \$1,370**

Code: IKAA25006458

## Newly Designed



## Even more reliability

It has a powerful and maintenance-free EC motor.

The surface is made of tempered and shatterproof safety glass, making it resistant to chemicals and easy to clean.

The housing is made of fireproof and robust die-cast aluminum with high quality powder coating.

Thanks to its enclosed design (protection class IP42), it can also be used under extreme conditions in the laboratory, for example when stirring cold samples.

The electronic connections on the back of the device are protected by the special design.

## Even more safety

Separate sensor for the safety temperature

Adjustable independent safety temperature

Hot surface warning

Confirmation mode

## For the budget conscious...

Same great IKA quality, the first RCT basic magnetic stirrer still offers reliability and safety. A long lasting partner in your lab.

**3 Units Only**

**Only \$690**

Code: IKAA20003581



Ref. Code:  
IKA-RCT

VALID UNTIL:  
31 Oct 2024  
for current stock only



1800 00 84 53

sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.