

**Miele**

Professional PG8504 Washer



Only \$7,833

Usually \$8,704

Code: MIEL62850422



Remove the monotony & hazard potential
of washing glassware manually,
let's put some precision and
safety into this procedure!

By removing the manual washing step from your wash
process there is reduced likelihood of glass breakage that
can result in cuts to lab staff.

- Powerful 500L/min pump with integrated heating elements enabling safe washing cycles with temperatures up to 70°C
- Underbench/Free-Standing
- Door dispensers for powder detergent and liquid neutraliser
- Suitable for operation on 240V 50Hz 13 amp supply
- 6 standard programs + 2 vacant program positions to meet your wash requirements
- Touch on steel controller with 3-line display for easy viewing of washer status
- Compact: 600mm x 600mm x 850mm W x D x H
- Stainless steel interior/exterior for durability and cleanliness

Miele - good for occupational health and safety!

Many options and accessories to select from.
Contact Bio-Strategy for assistance in configuring your washer.



Ref. Code:
MIEL-PG8504-1024

VALID UNTIL:
30 Nov. 2024



 1800 00 84 53

 sales.au@bio-strategy.com

 www.bio-strategy.com

 shop.bio-strategy.com