



# Vortex Mixer

Robust, Reliable, and Great Value!

## The Perfect Vortex Mixer For Your Everyday Mixing

**1 For \$450/each**  
**2 For \$425/each**  
**3+ For \$395/each**  
**Usually \$600/each**  
Code: OHAU30392118



The Ohaus Analogue Vortex Mixers offer variable speed control from low RPM start-up for gentle shaking to high-speed mixing for vigorous vortexing of samples.

Touch mode operation activates mixing when depressing the cup head or continuous mode for use with accessory attachments. The sturdy, rugged design ensures the unit stays in place while in use.

Each vortex mixer comes with a cup head and 3-inch accessory platform, to cater for a variety of vessels/tubes.

<b>Control:</b>	Analogue
<b>Orbit:</b>	4.9mm
<b>Speed Range:</b>	300 rpm – 2800 rpm
<b>Capacity:</b>	0.5kg
<b>Dimensions:</b>	8.3 in x 6.5 in x 4.8 in (LxHxW)
<b>Net Weight:</b>	4.4kg
<b>Warranty:</b>	2 years

### Applications

Re-Suspending Liquids, Mixing Reagents, Buffers and Emulsions, ELISA Assay, Drug Extractions, Mixing Samples

### Display

Independent LED displays for speed/time allow the operator to view both settings simultaneously (digital models) and adjustment knob with dial markings for variable speed control (analog models)

### Operation

Detachable 3-wire cord and plug (included)

### Communication

N/A

### Construction

Cup head and 3-inch accessory platform for beakers, flasks, test tubes etc. (included with all models), stainless steel microtube holder for 12 tubes (pulsing models)

### Design Features

Microprocessor-control to maintain set speed for consistent mixing action (digital models)



Part of DKSH Group

Ref. Code:  
OHAU-VM-1024

VALID UNTIL:  
31 Oct. 2024



1800 00 84 53

sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com