

Toploading Autoclave SX-700E 79 Litre

User Friendly | Ergonomic Height Design | Easy Operation



Only \$10,975

Was \$17,342

Code: TOMYSX-700E

Includes 2 stainless steel mesh baskets

**No plumbing necessary
Unit can be moved to where it is required!**

- Basket Dimensions: 345mm Ø x 181mm H
- Powered from a 240V 50Hz 15amp supply
- Simple lid opening using either hand or foot
- Grade 304 stainless steel for easy cleaning and durability
- Variable speed exhaust
- Three operating temperature ranges: Warming 45 to 95°C, Heating 45 to 104°C & Sterilizing 105 to 135°C
- Safety devices include: water level sensor, current leakage breaker, lid interlock, overheat prevention, over pressure prevention, open temperature sensor detection & safety valve
- Digital temperature & time display with analogue gauge for pressure
- Five sterilization profiles to select from (in addition to normal sterilization)
- Cooling fan to reduce run time
- Integrated (removable) condensate bottle to trap condensate from chamber during
- 2-year warranty

Good range of accessories also available, including punched metal stainless steel baskets, solid stainless steel baskets, hybrid basket: solid lower third with punched metal stainless steel upper, Durham tube stand, printer, external sensor (wandering probe) and preventative maintenance contracts.