



Wide Mouth Heavy Duty PP Bottle

Translucent, rigid polypropylene bottles with polypropylene screw caps, ideal for a wide range of liquid storage, media preparation and sampling applications

AZLOBWP0250P	AZLOBWP0500P	AZLOBWP1000P
250ml (10 items)	500ml (10 items)	1000ml (5 items)
Only \$16	Only \$18	Only \$19
Was \$24	Was \$28	Was \$30

- Leak-proof linerless cap design
- Smooth inner contours of neck for efficient pouring and ease of cleaning
- Suitable for food contact applications in accordance with EU and U.S. FDA regulations
- Good chemical resistance
- Autoclavable at 121°C
(NOTE: When autoclaving bottles always ensure the caps are completely loosened or removed to prevent accidental implosion)

Low Form Tapered Printed Graduations PP Beaker

Clear, blue, and easily readable graduations, along with a tapered design and ribbed supports, enable efficient stacking for storage optimization

AZLOBPM0250P	AZLOBPM0400P	AZLOBPM1000P	AZLOBPM2000P
250ml (10 items)	400ml (10 items)	1000ml (5 items)	2000ml (1 item)
Only \$12	Only \$17	Only \$13	Only \$4
Was \$20	Was \$30	Was \$23	Was \$7



- Manufactured from high clarity polypropylene
- Autoclavable at 121°C with good chemical resistance
- Precise non drip pouring spout supports efficient pouring
- *NOTE: Cleaning temperatures below 60°C are recommended to preserve marks and inscriptions*