

## 100 Litre Top-Loading Autoclave

**Only \$14,675**

Usually \$20,965

Code: TOMYFLS-1000

Includes two stainless steel mesh baskets & 2-year warranty

Powered from a 240V 50Hz 15amp supply and requires no plumbing so can be moved where required!

Large Capacity  
**100L**


Wide Chamber  
**45cm**



45cm Wide Chamber



- Large 100L vertical chamber capacity
  - Max Bottle Capacity: 34 (17 x 500mL per basket)
  - Max Flask Capacity: 4 x 5,000mL
- Integrated condensate bottle , no plumbing required!
- Full safety features including water level detection, and lid interlock
- Chamber dimensions: 450 Ø x 678 mm H
- External dimensions: 680 (D) x 680 (W) x 1022 (H) (including the pedal and control panel)
- Maximum sterilising temperature to 123 °C
- Weight (kg): 108



Easy opening foot  
pedal operated lid

Good range of accessories available including:  
punched metal stainless steel baskets, solid stainless steel baskets,  
printer, external sensor (wandering probe) and preventative maintenance contracts.