



PG8593 Under Counter Washer/Disinfector

Lab washer with liquid dispensing & DryPlus hot air drying




Only \$16,609

Usually \$18,054

Code: MIEL62858306

Included **FREE**
Second automatic dosing pump
for liquids plus 5L detergent and neutraliser
(valued at \$1,110)

- Single dosing pump for liquid detergent or liquid neutraliser, optional dispensing pump (free with this promotion) required for liquid neutraliser.
- Stainless steel interior/exterior for durability and cleanliness
- Spray arm and spray pressure monitoring ensuring a clean wash
- Compact: 600mm x 600mm x 850mm (820mm with lid removed)
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- Two wash levels
- Many options and accessories to select from
- Dispensing pump for liquid detergent,
- Programmable, 15 programs
- Demineralized water connection
- Capacity 128 flasks or 98 pipettes
- Multi-stage filtration system
- Intelligent A0-value controls
- 3 phase power, convertible to single phase (charges apply)
- Tested for 15,000 operating hours
- 2-year warranty



DryPlus

For perfect drying results

Active DryPlus HEPA filtered drying ensures your loads come out dry

Plus 15% discount on optional accessories – please call us to discuss your complete washing solution

Ref. Code: MIEL-PG8593

Valid Until: 31 Mar 2025

 **bio-strategy**
Part of DKSH Group



1800 00 84 53
sales.au@bio-strategy.dksh.com
www.bio-strategy.com
shop.bio-strategy.com

 **DKSH**

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.