



Adventurer™ Analytical Balance AX224

Only \$2,595

Usually \$3,591

Code: OHAU30100602

Includes a 3-year warranty

Intuitive balances designed for routine weighing, with a full-colour touch screen, automatic internal calibration (AutoCal), space-saving draft shield design, and 10 built-in applications.



Maximum Capacity: 220 g

Readability: 0.1 mg

Pan Size: 90 mm

Applications

Weighing, Parts Counting, Percent Weighing, Check Weighing, Dynamic Weighing, Formulation, Density Determination, Totalization, Display Hold, Batch Printing

Construction

Metal base, ABS top housing, stainless steel pan, glass draftshield with 2-piece top mounted side doors & sliding top door (1 mg models only), illuminated up-front level indicator, weigh below hook, security bracket, calibration lock, full housing in-use cover

Design Features

Selectable environmental filters, auto tare, user selectable span calibration points, software lockout and reset menu, user selectable communication settings and data print options, user definable project and user IDs, up to 9 operating languages



Protect your balance
with a dust cover.
Contact us for pricing.

Ref. Code: OHAU-AX224

Valid Until: 31 Mar 2025



1800 00 84 53
sales.au@bio-strategy.dksh.com
www.bio-strategy.com
shop.bio-strategy.com

