

Scientific Solutions

Trade In, Trade Up Your Old Viscometer

SPECIAL OFFER: UP TO 10% OFF



How old is your Brookfield viscometer?

Has it reached the end of life support from Brookfield?

DKSH and Brookfield are offering a discount up to 10% on new viscometer instrument purchases with the trade in of an old rotational viscometer. It can be any brand and does not need to be working!

Benefits of Trading In Your Old Instruments:

- Upgraded Technology
- Accuracy and Precision
- Efficiency and Productivity
- Reduce Downtime
- Cost Saving
- Data Management
- Competitive Advantage

Product Family	End of Repair/ Service (EOS)	End of Life (EOL)	Recommended Alternative
Dial/Analogue			DVPlus
DV-E			DVPlus
DV-I			DVPlus
DV-I*			DVPlus
DV-I Prime			DV2T or DVNext
DV-II			DV2T or DVNext
DV-II*	July 1, 2020	Dec. 31,2022	DV2T or DVNext
DV-II Pro			DV2T or DVNext
DV-II* Pro			DV2T or DVNext
DV-II Pro Extra			DV2T or DVNext
DV-III			DVNext
DV-III*			DVNext
DV-III Ultra			DVNext

* Discount applies to instruments only not accessories or software,
 * One trade in instrument per new instrument purchase.
 * Please include Trade In on your PO to receive this discount.
 * This offer cannot be combined with existing discounts or promotions.

Send in Your Trade-In Request Today!

Offer Valid Until:
31 Aug 2025



Featured Viscometers



DVPlus Viscometer

Features:

- 5-inch Full-Color Touch Screen Display
- Advanced interface with enhanced controls and a real-time test indicator
- USB and Bluetooth* interface for use with the DVPlus Connect app
- Temperature offset capability to $\pm 5^{\circ}\text{C}$
- Choice of multiple languages, including English, French, German, Portuguese, Chinese, Japanese, and Spanish
- Total of 740 speeds for an excellent range capacity
- Optional RTD Temperature Probe (DVP-94Y)
- Direct access to the time measurement function (time to torque, time to temperature, time to stop).
- Automatic range calculation
- Accuracy of $\pm 1.0\%$ of range
- Repeatability of $\pm 0.2\%$
- 17025 Calibration (Optional)
- Model G Base



DVNext Rheometer

Features:

- Digital Leveling
- Start-up Wizard
- Gel-Timer Functionality
- Optional Magnetic Coupling Attachment (replaces EZ-lock)
- Optional Bar Code Scanner
- 17025 Calibration (Optional)
- Hand-held Scanner for Bar Code Detection
- Accessory Detection with Bar Code Detection
- Magnetic Coupling Attachment (replaces EZ-Lock)
- Compliance to 21CFR Part 11 in Stand-alone and GAMP Compliant
- Ethernet Connectivity
- LIMS Connectivity
- Printing to non-editable PDF



DV2T Touch Screen Viscometer

Features:

- 5-inch full color, touch screen display supports multiple languages
- Displayed info includes: viscosity (cP/mP·s), temperature ($^{\circ}\text{C}/^{\circ}\text{F}$), shear rate/stress, % torque, spindle/speed, step program status
- Enhanced security provides customizable user levels and password access.
- Built-In options include timed tests, data averaging, programmable QC limits/alarms, customizable speed/spindle lists, on screen data comparison
- Auto range shows maximum viscosity measured with any spindle/speed combination
- USB PC interface provides optional computer control and automatic data gathering capability Download custom test programs with included PG Flash software
- Front-facing bubble level for convenient viewing
- Built-in temperature probe
- Accuracy of $\pm 1.0\%$ of range with displayed test data
- Repeatability of $\pm 0.2\%$
- 17025 Calibration (Optional)
- NIST traceable viscosity standards available
- Also available: Magnetic coupling system, durable ball bearing suspension system, quick action lab stand and RheocalcT software.

Offer Valid Until:
31 Aug 2025

For more information on Brookfield Viscometers, or to discuss your requirements, contact us now:

Phone: 1800 010 114

Email: customerservice.ins@dksh.com

Web: labshop.dksh.com.au/shop/



Visit DKSH LabShop Now!